
SINTERFILTER
VitraPOR

ROBU[®]
GLASFILTER-GERÄTE

TECHNISCHE DATEN
TECHNICAL DATA

Der Name

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The Name

VitraPOR is the trade-name for a complete series of glass filter-products, equipped with filter elements. ROBU-products represent the internationally established type of borosilicate glass 3.3 according to DIN-ISO 3585 and english, american and french standards. Highest possible chemical resistance, minimum thermal expansion as well as the high thermal shock resistance are the characteristic points of ROBU glass filter products. Due to these features the ROBU filters are specially designed for use in the chemical, biotechnical, pharmaceutical and laboratory field.

Standard fritted discs are available:

- in round shape from 5 mm to 120 mm diameter
- in round shape bi-plane ground
- in round shape bi-plane ground and with fused edge
- filter candles in cylindrical or conical shapes from 9 mm to 40 mm diameter
- all filterapparatus and absorbers

Custom made **VitraPOR**-filters:

- any shape up to 400 mm in almost any thickness
- all custom made discs are available unground, bi-plane ground and with fused edge
- any cylindrical or conical filter candle according to your design

Please contact us for any custom made item. We are the experts for any special design and are able to solve your problems fast and at reasonable costs.

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
ROBU VitraPOR filter

The filtering effect of a porous body consists of the sieve effect and the adsorption. If a filter is to serve as a sieve only, the filter-pores have to be smaller than the smallest solid particles to be filtered. Therefore the amount of liquid passing through is relatively small. Larger amounts of liquid can be achieved with a filter body, whose filtering effect is based upon adsorption. In this case the filter pores can be considerably larger than the solid particles to be separated. During sterile filtration of air, for example, with filter bodies of a pore-size of 15-40 microns, the size of the particles to be retained is only a 1/10 to 1/50 of the pore-size. A safe adsorption, however, can only be achieved when the inner surface of the filter-body is very large. In order to obtain this large inner surface, any single grain has to have a maximum surface area. **VitraPOR** filters are produced in a sophisticated manufacturing process. Special milling procedures are used to obtain a splittling grain with an extra-large surface. With a certain post-production treatment this surface is additionally enlarged.

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Introduction

Borosilicate glass 3.3 is specified as a standard glass - it contains mainly sand, calcium carbonate and sodium carbonate.

High quality sand, dehydrated borate boric acid, alumina and salt are used in its manufacture and strict control is applied to ensure the specifications for purity and consistency of these materials are maintained. Because of its resistance to heat, corrosion and thermal shock borosilicate glass 3.3 is used extensively in science and in industry.

Chemical Composition

The following is a typical analysis of borosilicate glass 3.3:

Element	% by weight
Silica (SiO ₂)	80.60
Boric oxide (B ₂ O ₃)	12.60
Sodium oxide (Na ₂ O)	4.20
Alumina (Al ₂ O ₃)	2.20
Iron oxide (Fe ₂ O ₃)	0.40
Calcium oxide (CaO)	0.10
Magnesium oxide (MgO)	0.05
Chlorine (C1)	0.10

There may be heavy metals in a concentration of less than 5 ppm.

Physical Properties

Coefficient of Expansion	33 x 10 ⁻⁷ /°C	between 20°C - 300°C
Specific Heat	750	J/KG°C bei 20°C
Thermal Conductivity	1.13	W/m²K at 20°C
Density	2.23 x 10 ³	kg/m ³
Poisson's Ratio	0.22	between 25°C - 400°C
Young's Modulus	6.500	kg/mm ² at 25°C
Rigidity Modulus E	62,5	kN/mm ²
DPH (Vickers) Hardness	580	Kg/mm ² with 50 gram load
Relative Hardness	1.52	(comparative Soda-Lime = 1.0)
Refractive Index	1.474	Sodium D - Line
Dielectric Constant	4.6	at 1 MHz and 20°C
Loss Factor	2.6	% at 1 MHz and 20°C
Log ₁₀ Volume Resistivity	15	Ohm - cms at 20°C
Surface Resistivity	10 ¹³	Ohms/cm ² at 50% humidity

Thermal Expansion

The coefficient of linear thermal expansion is one of the characteristic properties by which glass is measured and is defined as the change in unit length per degree rise in temperature. The thermal expansion is 33 x 10⁻⁷ per °C between 20°C and 300°C and is controlled to a tolerance of ± 0.5 x 10⁻⁷ per °C. Due to its porous structure, sintered glass has a slightly different thermal expansion. Care should be taken when fusing it to massive borosilicate glass 3.3.

Thermal Resistance

Since borosilicate glass has a high softening point, it can be used at temperatures of the order of 515°C. Care must be exercised, however, in cooling down apparatus which has been held at high temperatures. The cooling process must be slow and uniform, particularly in the case of thick sections (see annealing). If a temperature gradient is applied to the glass, stresses are set up - the hotter glass being under compression and the cooler under tension. The magnitude of the tensile stress depends on the temperature difference which in turn depends on the glass thickness. ATTENTION: Abrasions greatly reduce the temperature shock resistance.

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
Maximum Temperatures

In general it is recommended, that the strain point (515°C) be regarded as the maximum safe operating temperature of borosilicate glassware 3.3. For some shapes and for a relatively time this limit can be exceeded, but at 580°C there is danger of deformation and in case of sintered glassware the porous structure may be changed. At high temperatures the glass may acquire permanent stress on cooling and this may result in subsequent breakage.

If it is suspected that permanent stress has occurred, the article should be annealed making references to the annealing process explained below. Permanent stress can greatly reduce the mechanical and thermal resistance.

Viscosity

Although glass has no defined softening point, there are four points in the viscosity-/temperature relation which have accepted definitions:

1. Strain Point 515°C
Viscosity of the glass is $10^{13.5}$ Ns/m² ($10^{14.5}$ dPas)
2. Annealing Point 565°C
Viscosity of the glass is $10^{12.4}$ Ns/m² ($10^{13.4}$ dPas)
3. Softening Point 820°C
Viscosity of the glass is $10^{6.6}$ Ns/m² ($10^{7.6}$ dPas)
4. Working Point 1225°C
Viscosity of the glass is $10^{3.0}$ Ns/m² ($10^{4.0}$ dPas)

Annealing

The annealing of glass is the process by which it is heated to and held at a controlled temperature for a defined period to relieve internal stresses. Careful cooling under controlled conditions is essential to ensure, that no stresses are reintroduced by chilling.

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Tabelle / Table

Refractive index (sodium D line)	1.474
Dispersion (NF-NC)	0.00738
Visible light transmission	92% in 2 mm thick glass
.	91% in 6 mm thick glass
Stress optical coefficient	3.5 nm/kg cm ² per cm length

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at 20 °C 1.13 W/m°K
 150 °C 1.25 W/m°K
 300 °C 1.42 W/m°K

Chemical Resistanc

The chemical resistance of glass is better than that of all other known materials. It shows resistance to water, acids and alkalis, salt, organic substances as well as chlorine and bromine. Hydrofluoric acid, concentrated phosphoric acid and strong alkalis attack the glass surface at higher concentrations and temperatures.

The following are typical results of tests undertaken to international standards:

1. Hydrolytic resistance
 Glass grains (300 - 500 microns)
 acc. ISO 719 (98°C) Water attack: < 0,01 mg Na₂O/g grains
 Class ISO 719 - HGB 1
 acc. ISO 720 (121°C)
 Class ISO 720 - HGA 1

2. Acid Resistance
 Resistance acc. DIN 12116 Acid Class 1
 Resistance acc. DIN/ISO 1766 = Loss in weight < 100 µg Na₂O/dm²

3. Alkali Resistance
 Resistance acc. DIN 52322 / ISO 695 Class A2

Optical Properties

Massive borosilicate glass 3.3 has a high ability to transmit light throughout the visible range of the spectrum. Ranging from approx. 310 to 2300 nm there is virtually no absorption through the thickness of glass generally encountered in laboratory ware (up to approx. 6 mm). The following table shows additional data:

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






Porosity

The VitraPOR glass filters are divided into seven porosity classes from 1,0 µm to 500 µm. Therefore they are suitable for all applications in the chemical and pharmaceutical industry, in laboratories as well as in large technical plants.

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
Druckunterschied/Pressure differential Δp (mbar)

Flow rate and Pressure Drop

Flow rates of liquids and drop in pressure of gases are functions of:

- diameter of the filter
- porosity of the filter
- relation of pressure/vacuum in the apparatus
- physical properties of the medium to be filtered.

To determine the possible application of glass filter-products and to choose the best apparatus for the work intended, it is necessary to know the flow rates of liquids and gases. Please use our diagrams 1 and 2 which are shown for water and air. The data apply to filter discs of 30 mm diameter with average properties and usually effective filter area and - thicknesses, as well as poresizes.

The flow rates for other diameters can easily be taken by multiplying this factor by the conversion factor. Minor deviations with very limited tolerances are possible. Strict control is recommended, especially when filtering bacteria using porosity 5.

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Cleaning

Successful experimental results can only be achieved by using clean apparatus. The following advice may be used for the massive and porous borosilicate glass 3.3:

Wash glassware as quickly as possible after use. During washing with water all parts of the article should be scrubbed with a soft brush. Avoid any abrasion of the glassware! After cleaning thoroughly rinse with distilled water. Drying can be undertaken either in air or at a temperature not exceeding 100°C.

If after mechanical cleaning for example filter-pores should still remain clogged, the thorough chemical cleaning is required. A few suggestions found generally useful follow:

Contamination

Solvent

Barium sulphate

Hot concentrated sulphuric acid (100°C) Silver-chloride Ammonian or Sodium hyposulfite

Copper or Iron Oxides

Hot Nitric acid

Mercury residue

Hot Hydrochloric acid plus potassium chlorate

Mercury sulphide

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Albumen

Hot hydrochloric acid or hot ammonia

Grease, oil, fatty materials

Carbon tetrachloride

Organic matters

Hot concentrated cleaning solution, or hot concentrated sulfuric acid plus a few drops of sodium or potassium nitrite

Glucose

Carefully heat to approx. 200°C with mixed acid (5% H₂SO₄ and 1% HNO₃)

Extensive rinsing with water must be obviously follow.

Warning!

Hot concentrated phosphoric acid and hot alkaline solutions attack the glass surface, they are unsuitable as cleaning agents. If they have to be filtered, an increase in pore size and reduced life of the fritted disc is unavoidable.

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Personal Safety

When treated with proper care borosilicate glass 3.3 laboratory apparatus will give long and satisfactory service. It is ecologically harmless and can be reused and recycled many times. The following notes have been prepared to assist users in obtaining the maximum life from their apparatus:

Initial Cleaning

Before the first use of glass-filters, they should be cleaned with diluted, warm hydrochloric acid, followed by several rinses with distilled water. Glassfilters should be cleaned immediately after use. The surfaces can be wiped with a squeegee or brush. Surface rinsing with a wash bottle or backflushing under water tap is often sufficient.

Pressure and Vacuum

Even small scratches or internal abrasions of the glass surface can greatly reduce the vacuum- or pressure resistance. Filterpores may clog during filtration and cause a pressure drop - in any case 1 kp/cm² may not be exceeded. Never subject glassware to sudden pressure changes.

Temperature Resistance

The maximum allowable service temperature is approx. 500°C. To avoid permanent stress, careful heating and cooling according to the schedule for borosilicate glass 3.3 should be ensured. Wet glass-filters should slowly be heated up to 80°C and dried for one hour, before temperature is increased. Permanent stress can reduce the mechanical and thermal resistance. Never subject glassware to sudden temperature changes.

Chemical Resistance

High concentrations of hydrofluoric acid or hot alkali solutions attack the glass and can cause a deterioration in the filter and increase of the pore size.

Safety Advice

When working with glassware always wear protective glasses and gloves to avoid injuries. Use safety screen, -hood or similar protection when working with glassware subject to pressure or vacuum. Be prepared, that reagents may leak from a broken vessel.

Please heed these guidelines in combination with the respective statespecific regulations for the use of glassware in the laboratory.

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The Twin-Filter...

...combines two different porosities in a single filter, made from pure borosilicate glass 3.3. A coarser layer serves as a pre-filter and prevents the filter pores from clogging.

The finer part holds back even smaller particles and guarantees an optimal filtration result.

The built-in pre-filtration feature extremely extends the life of the filter. In continuous processes, the performance will be considerably better, while the fine pores are well protected by the mechanical strength of the coarser side.

The Twin-Filters are made by the proprietary ROBU high temperature Hard-Sintering process. As with all ROBU filters, no binders or additives are allowed in any step of the manufacturing process. This is your assurance, that the Twin-Filter meets the highest laboratory standards for purity. Different combinations of pore sizes acc. ISO 4793 are available to meet the requirements in chemical, pharmaceutical, biological and medical processes.

Upon request, we are able to produce special designs and other porosity-combinations to meet your unique requirements.

The Twin-Filters are available as blanks in round or rectangular shape from 5 mm to 400 mm diameter, or fused into all standard filter-apparatus such as funnels, crucibles and bubblers.

Porosity	Pore Size	
	Coarse	Fine
0/2	160-250 µ	40-100 µ
1/3	100-160 µ	16- 40 µ
2/4	40-100 µ	10- 16 µ
3/5	16- 40 µ	1- 1,6 µ

Special designs

Based upon the wide-spread pallet of standard glassfilters, ROBU was able to build up a good reputation on the international market. Special designs glass filters and -apparatus are available in almost any shape, size and porosity for many extraordinary applications.

The high quality-standards of our production, the engagement in developments and the interaction with our customers show our potential for the chemical, pharmaceutical, medical and environment industry.

If you require further information on ROBU glass filters, their specifications and applications, please don't hesitate to contact our qualified staff.

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Our wide range of porosities and sizes makes it easy for you. If your specific demands cannot be met with our standard products, we will work closely with you to create what you need. ROBU can fabricate sinterfilters in almost any shape and size.

ROBU®

ROBU GLASFILTER-GERÄTE GMBH
Schützenstr. 13, D-57644 Hattert, Germany
Tel.: 02662/8004-0
Fax.: 02662/8004-40
info@robuglas.com
www.robuglas.com